


# KWIK KLEEN<sup>®</sup>

## KERRISON

### FEATURES:

- Patented design enables the Kerrison to be fully cleaned, thus eliminating the build up of bone, tissue, and bio-burden that can cause sticking and non-performance.
- Wide grip handles provide complete comfort and control.
- Single lever makes ease of opening and closing.
- Bite size numbers are etched on both sides of Kerrisons.
- Proven design provides a consistent smooth cutting action.
- 180 mm - 200 mm - 230 mm length shaft available - also custom lengths available.
- Top rail does not detach, thus eliminating errors in assembly.
- Manufactured from 100% high grade German stainless steel.



Advanced  
Endoscopy  
Devices

*Endoscopy Specialists since 1985*

22134 Sherman Way  
Canoga Park, CA 91303 - U.S.A.  
phone: 818-227-2720 - fax: 818-227-2724  
www.aed.md  
E-mail: advendodev@earthlink.net

# KWIK KLEEN<sup>®</sup>



## KERRISON

MODEL #	DESCRIPTION:		
	Bite Size	Length	Angle of Tip
8527-00	1 mm	x 180 mm	x 40 Degree
8527-01	2 mm	x 180 mm	x 40 Degree
8527-02	3 mm	x 180 mm	x 40 Degree
8527-03	4 mm	x 180 mm	x 40 Degree
8527-04	5 mm	x 180 mm	x 40 Degree
8527-A6	6 mm	x 180 mm	x 40 Degree
8527-05	1 mm	x 200 mm	x 40 Degree
8527-06	2 mm	x 200 mm	x 40 Degree
8527-07	3 mm	x 200 mm	x 40 Degree
8527-08	4 mm	x 200 mm	x 40 Degree
8527-09	5 mm	x 200 mm	x 40 Degree
8527-B6	6 mm	x 200 mm	x 40 Degree
8527-10	1 mm	x 180 mm	x 90 Degree
8527-11	2 mm	x 180 mm	x 90 Degree
8527-12	3 mm	x 180 mm	x 90 Degree
8527-13	4 mm	x 180 mm	x 90 Degree
8527-14	5 mm	x 180 mm	x 90 Degree
8527-C6	6 mm	x 180 mm	x 90 Degree
8527-15	1 mm	x 200 mm	x 90 Degree
8527-16	2 mm	x 200 mm	x 90 Degree
8527-17	3 mm	x 200 mm	x 90 Degree
8527-18	4 mm	x 200 mm	x 90 Degree
8527-19	5 mm	x 200 mm	x 90 Degree
8527-E6	6 mm	x 200 mm	x 90 Degree



**Advanced  
Endoscopy  
Devices**

*Endoscopy Specialists since 1985*

22134 Sherman Way  
Canoga Park, CA 91303 - U.S.A.  
phone: 818-227-2720 - fax: 818-227-2724  
www.aed.md  
E-mail: advdodev@earthlink.net