




**Advanced  
Endoscopy  
Devices**

# LED Stroboscope



**8000-S**

- Vocaleled 5.0
- LED Compact Stroboscopic Light Source
- Ideal for laryngostroboscopic practices

**LED Handpiece**

- White LED Lamp, Comparable 100 Watt
- 6500 K (+/- 500K)
- Lifetime: 20000 hr
- Connector Type: Storz

**Features**

- Stroboscope light / continuous light
- Stroboscope range: 30 - 1200 Hz
- Automatic stroboscope function
- Slow motion rate: 0,5 / 1 / 1,5 / 2 Hz
- Brightness control




**LED Handpiece**



**Stethoscope**



**Foot pedal (optional)**



- Excellent Color Contrast
- High Resolution
- Wide Angle

**8780C**

- 70° HD Laryngoscope
- 8mm x 18cm

**8780D**

- 90° HD Laryngoscope
- 8mm x 18cm

**1-818-227-2720**

**Advanced Endoscopy Devices, Inc.**  
22134 Sherman Way, Canoga Park, CA 91303 USA  
www.aed.md sales@aed.md