



Advanced  
Endoscopy  
Devices



## FLEXIBLE ENDOSCOPE FEATURES

- Naso-Pharyngo-Laryngo-Fiberscope
- TRU - VU® high resolution fiber optic imaging system
- Two Models: 2.8mm and 3.4mm
- 0 degree Angle of View
- 150 degree distal deflection (both ways)
- Large 80 degree Field of View
- Fully immersible
- Gas Sterilization
- Depth of Field 1 – 50mm
- Includes carrying case, ETO-Cap, leak tester and instruction manual

### Model 801.3400

0° AOV  
3.4mm x 30cm  
150° Deflection (both ways)

### Model 801.2801

0° AOV  
2.8mm x 30cm  
150° Deflection (both ways)



**1-818-227-2720**

Advanced Endoscopy Devices, Inc.  
22134 Sherman Way, Canoga Park, CA 91303 USA  
[www.aed.md](http://www.aed.md) [sales@aed.md](mailto:sales@aed.md)