



Advanced  
Endoscopy  
Devices

# Bipolar Saline Resectoscopes

## AED Bipolar Working Elements 24/26 fr

- TUR in Saline (NaCl 0.9%)
- No Need For Patient Plate
- Precise Cut and Coag.
- Available with Passive or Active Working Element
- Wide Array of Electrodes Available
- 24/26fr Slim CFR Sheath
- Light and Well Balanced Design
- AED Kwik Snap Sheath Mount


32-4220BW  
Passive



32-4216BW  
Active



## AED Saline Resectoscope 24/26 fr



1-818-227-2720

Advanced Endoscopy Devices, Inc.  
22134 Sherman Way, Canoga Park, CA 91303 USA  
[www.aed.md](http://www.aed.md) [sales@aed.md](mailto:sales@aed.md)