



**Advanced
Endoscopy
Devices**

Tru-Vu Hysteroscopes

Diagnostic



5729BL

- 2.9mm Diameter
- 30 Degree
- Hi-Resolution
- Excellent color reproduction



8805B

- 4mm Diameter
- 30 Degree
- Hi-Resolution
- Excellent color reproduction



HD2304

- Hystero Diagnostic Sheath
- 30 Degree
- For 2.9mm Scope
- 4mm Diameter



HD2305

- Hystero Diagnostic Sheath
- 30 Degree
- For 4mm Scope
- 5mm Diameter



HD2304-CF

- Hystero Diagnostic Sheath
- 30 Degree
- Continuous Flow
- For 2.9mm Scope
- 5mm Diameter



HD2305-CF

- Hystero Diagnostic Sheath
- 30 Degree
- Continuous Flow
- For 4mm Scope
- 6.5mm Diameter

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md