


**Advanced
Endoscopy
Devices**

Mini-Op Laparoscopy

Semi-Rigid Forceps




308-013-0385H

- Alligator Fenestrated Grasper
- Double Action
- 5 FR Diameter
- 38cm Working Length



300-012-0380H

- Biopsy Forcep
- Double Action
- 5 FR Diameter
- 38cm Working Length



300-016-0380H

- Blunt Scissor
- Single Action
- 5 FR Diameter
- 38cm Working Length



308-013-0385H2

- Smooth Fenestrated Grasper
- Double Action
- 5 FR Diameter
- 38cm Working Length



300-017-0380H

- Hook Scissor
- Single Action
- 5 FR Diameter
- 38cm Working Length



294-016-0380H

- Pointed Scissor
- Single Action
- 5 FR Diameter
- 38cm Working Length

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md