



Advanced
Endoscopy
Devices

SHAVER *Leak Tester*



Designed for pressure testing arthroscopic shavers to help identify leaks caused by failing seals that degrade over time from repeated use and processing.

The Shaver Leak Tester assists healthcare workers in reprocessing areas to make sure their shaver is working properly by using both a hand pump to apply pressure and a testing stop made from polypropylene to create a temporarily closed system within the fluid pathway of the shaver to pressure test for leaks.

Reduce the potential for cross-contamination, damage and costly repairs that result from using a shaver with leaks

Available for the following shavers:

- Arthrex shaver
- Conmed Linvatec shaver
- Dyonics shaver
- Midas Rex shaver
- Stryker Formula shaver
- Stryker Small Joint shaver



1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md