



ATLAS[®] PIN & CABLE CUTTERS


FEATURES:

- Double the Power of a Standard Pin Cutter
- Unique Leverage Force Torque System
- 100 % German Stainless Steel Construction


Model # 6500
Extractor Plier Forcep




Model # 6300
Pin - Wire
End & Side Cutter
max. 1.6 mm



Model #6400
Cable Cutter
for cables up to 2.0 mm



Model #6200
Pin - Wire Side Cutter
max. 3.0 mm
with Rubber Jaw Pads



- Replaceable Tungsten Carbide Blades
- Replaceable Rubber Pads
- TiN Coated Jaws
- Steam Autoclavable



Advanced
Endoscopy
Devices

Endoscopy Specialists since 1985

22134 Sherman Way
Canoga Park, CA 91303 - U.S.A.
phone: 818-227-2720 - fax: 818-227-2724
www.aed.md
E-mail: advendodev@earthlink.net