



Advanced
Endoscopy
Devices

Precise™

Lap instruments for every budget

Precise line of instruments were designed with the simplicity of first generation laparoscopic instrumentation, yet with the innovations of third generation instruments. In other words they are sophisticated and easy to use. **Precise** lap forceps completely disassemble for thorough and effective cleaning. These instruments were made for **reprocessing ease**. The **Precise action** of the jaws offer excellent tactile feedback and the **large rotator** allows for **Precise** positioning.

- 1st generation simplicity
- 3rd generation innovation
- Large rotator
- Color coded rotator caps
- Affordable price
- Reprocessing ease
- Precise jaw action
- Excellent tactile feedback
- 45 popular jaw patterns
- 5mm and 10mm diameters
- 33cm and 45cm lengths
- Standard or ratchet handle options available
- Highest quality stainless steel jaws
- Sold and serviced in the USA
- FDA Registered & ISO13485 Certified

Jaws
Fully Disassemblable



Rotating Knob
Controlled 360° rotation



45° Caution Pin
Reduces electrocautery cord interference



Insulated Shaft
Arc resistant polymer



Optional
Color-Coded Lock Caps



Shaft-Handle
Easily disassembled



Economically, you won't find laparoscopic instrumentation that better fits the demand of conventional standards and your budget. With **four times the savings** compared to other laparoscopic instrumentation, and durability necessary for repetitious use, you'll find it hard to dispute the value of **Precise** instruments. **Precise** laparoscopic instrumentation are the instruments for every budget.

Order yours today!
1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md