



Advanced
Endoscopy
Devices

SPEEDLOCK[®]

Surgeon's Choice

AED's third generation **SPEEDLOCK** modular forceps are the industry's premier laparoscopy forceps. Its quality and handling surpasses any other product line in the endoscopic market today. With nearly **4,000** possible combinations, **SPEEDLOCK** comes in **4** shaft lengths, **2** diameters, more than **100** jaw patterns, and **10** handle options. With AED's **SPEEDLOCK** line of instruments there are **unlimited variations**.

- Modular Laparoscopic Instrument
- Colored ID rings for rotator
- Flushable
- Unique jaw bayonet coupling
- Over 100 popular jaw patterns
- 5mm and 10mm diameters
- 26/30/34/45cm lengths
- 10 popular ergonomic handle choices available
- High-grade German stainless steel jaws
- Environmentally friendly

Jaws
Features a unique bayonet locking system



Flush Port
Easily remove debris and reprocess



45° Cautery Pin
Reduces electrocautery cord interference



Insulated Shaft
Arc resistant polymer



Optional
Colored rings for rotator



Shaft-Handle
Easily disassembled



With our years of experience and continuous development, we have made this exceptional modular instrument line possible and that's why using AED's **SPEEDLOCK** instruments means being one step ahead.

Order yours today!
1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md