

Precise™









Lap instruments for every budget

















FEATURES








- *Modular 3 Piece Design*
- *Precise Jaw Action*
- *360° Rotation – Large Rotator*
- *Fully Disassembled For Ease of Cleaning*
- *Excellent Tactile Feedback*
- *Cost Effective – Competitively priced*
- *46 Popular Jaw Patterns*
- *Available in 5mm and 10mm Diam.*
- *Available in 33cm and 45cm Length*
- *Surgeon's Choice*
- ***1st Generation Simplicity 3rd Generation Innovation***

















	<p>P-1290</p>	<p>Retraction Forcep Double Action 5mm</p>
	<p>P-1295</p>	<p>Retraction Forcep Single Action 5mm</p>
	<p>P-1300</p>	<p>Blunt Grasping Forcep 5mm</p>
	<p>P-1301</p>	<p>Fenestrated Grasping Forcep 5mm</p>
	<p>P-1302</p>	<p>Maxi Grasper 5mm</p>
	<p>P-1303</p>	<p>Croc Jaw Grasper 5mm</p>






	<p>P-1312</p>	<p>Fenestrated Curved Grasper Single Action 5mm x 22mm Jaw</p>
	<p>P-1313</p>	<p>Fenestrated Curved Grasper Double Action 5mm x 32mm Jaw</p>
	<p>P-1317</p>	<p>Fine Dissector 5mm</p>
	<p>P-1318</p>	<p>Maryland Jaw 5mm</p>
	<p>P-1323</p>	<p>Endoclinch Forcep 5mm</p>
	<p>P-1336</p>	<p>Mixer Clamp 10mm</p>
	<p>P-1337</p>	<p>Allis Grasper 5mm</p>

	P-1340	Maryland Grasper 10mm
	P-1343	Bowel Grasper 5mm
	P-1354	Standard Grasper 5mm
	P-1375	Wave Grasper Double Action 5mm
	P-1376	Mixer Clamp Right Angle 5mm
	P-1377	Retention Forcep 5mm
	P-1378	Adhesion Grasper 5mm

	P-1379	Fenestrated Wave Grasper 5mm
	P-1380	Claw Grasper 5mm
	P-1381	Wave Grasper Single Action 5mm
	P-1401	Cobra Grasper 5mm
	P-1411	Kelly Forcep 10mm
	P-1413	Straight Dissecting Forcep 5mm
	P-1431	Hook Scissor 5mm

	P-1459	Straight Scissor Single Action 5mm
	P-1460	Straight Metz Scissor Double Action 5mm
	P-1461	Curved Metz Scissor Double Action 5mm
	P-1485	Biopsy Spoon Forcep 5mm
	P-1490	Biopsy With Tooth Forcep 5mm
	P-1536	Tenaculum 10mm
	P-1537	Tenaculum 5mm

	<p>P-1541</p>	<p>Babcock Forcep 5mm</p>
	<p>P-1545</p>	<p>Babcock Forcep 10mm</p>
	<p>P-1548</p>	<p>Short Wave Grasper 5mm</p>
	<p>P-1560</p>	<p>Bowel Grasper 10mm</p>
	<p>P-1580</p>	<p>Duval Grasper 5mm</p>
	<p>P-1605</p>	<p>Claw Forcep Single Action 5mm</p>
	<p>P-1606</p>	<p>Claw Forcep Double Action 10mm</p>

	<p>P-1607</p>	<p>Claw Forcep Single Action 10mm</p>
	<p>P-1660</p>	<p>Spoon Forcep 10mm</p>
	<p>P-8514</p>	<p>Johann Double Action 5mm x 22mm Jaw</p>
	<p>P-JRT-10</p>	<p>Hunter Grasper 5mm</p>
	<p>P-MDF-100</p>	<p>Medifex Grasper 5mm x 55mm Jaw</p>

	<p>P-1356</p>	<p>Shaft For Precise Forcep 5mm x 33cm</p>
	<p>P-1356-45</p>	<p>Shaft For Precise Forcep 5mm x 45cm</p>
	<p>P-1357</p>	<p>Shaft For Precise Forcep 10mm x 33cm</p>
	<p>P-1357-45</p>	<p>Shaft For Precise Forcep 10mm x 45cm</p>

 <p>A dark metal laparoscopic handle with a curved, hook-like tip and two circular finger loops. It has a small white dot on the upper handle and a threaded end at the top.</p>	<p>P-1383</p>	<p>Standard Precise Handle</p>
 <p>A dark metal laparoscopic handle with a curved, hook-like tip and two circular finger loops. It features a ratchet mechanism on the lower handle, indicated by a small white label. It has a small white dot on the upper handle and a threaded end at the top.</p>	<p>P-1384</p>	<p>Ratchet Precise Handle</p>

	<p>P-RC100</p>	<p>Red Rotator Cap for Precise Forcep</p>
	<p>P-RC200</p>	<p>Blue Rotator Cap for Precise Forcep</p>
	<p>P-RC300</p>	<p>Green Rotator Cap for Precise Forcep</p>
	<p>P-RC400</p>	<p>Yellow Rotator Cap for Precise Forcep</p>
	<p>P-RC500</p>	<p>Gray Rotator Cap for Precise Forcep</p>
	<p>P-RC600</p>	<p>Black Rotator Cap for Precise Forcep</p>