




Advanced  
Endoscopy  
Devices

# Meeker Dissector Forceps



Meeker Dissecting Forceps  
AED116850



- Monopolar
- 360° Rotation
- Flushable port
- The 90° curve provides excellent visualization
- 10mm x 33cm
- **Ratchet Available \*AED116800**

Mixer Dissecting Forceps  
AED116703



- Monopolar
- 360° Rotation
- Flushable port
- The 60° curve gives optimum visualization of duct and artery
- 10mm x 33cm
- **Ratchet Available \*AED116732**

**1-818-227-2720**

Advanced Endoscopy Devices, Inc.  
22134 Sherman Way, Canoga Park, CA 91303 USA  
www.aed.md sales@aed.md