



Advanced
Endoscopy
Devices

Tru-Vu Laser Resectoscope

Laser Resectoscope System

(For laser fibers up to 1.5mm or 1500µm)



Model # 32-4220LR

- Laser Resectoscope Working Element
- Passive Action System
- Ergonomic Design Handle
- Large 1500µm laser channel
- LRS works with any brand of fiber laser
- Modular system is compatible with existing AED or Storz style equipment
- Inner Sheath allows surgeons to change instruments without removing the outer sheath during a case

* **Scope & Sheath not included with Resectoscope**

Additional Items Available

32-4205RR - 26FR, Resectoscope Sheath

8805B - 30 Degree, AED Gold Cystoscope

150050L - Telescope Bridge, with 7FR Channel

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md