


Tru-Vu[®] Continuous Flow Laser Cystoscope 23 FR

FEATURES:

- 23 FR Small Diameter Sheath
- Large Port Continuous Flow System
- 7 FR Working Channel
- Compatible for Green Light[™] PVP Procedure*
- Profile Distal Tip Designed for Easy Passage
- Tru-Vu[®] Precision Rod Lens Cystoscope
- 100% Stainless Steel Construction
- Ease of Reprocessing - Cold Soak - ETO Gas - Steam Autoclavable

8805B- 4.0 mm x 30 Degree Tru-Vu[®] Rod Lens Cystoscope. Steam Autoclavable Color Code: Red



32-4085 LS - 23 FR Continuous Flow Laser Cystoscope Sheath with Hi Flow Ports

32-4085 LO - Obturator for 23 FR laser Cystoscope Sheath

32-4085 LB - Continuous Flow Laser Cysto Bridge and Inner Sheath with 7 FR Working Channel

32-4085 VO - Visual Obturator for 23 FR Laser Cysto Sheath



32-4120 - Scissor Forcep 7 FR Flexible




32-4130 - Grasping Forcep 7 FR Flexible



32-4135 - Biopsy Forcep 7 FR Flexible

*Green Light[™] is a trademark of AMS Inc.



Advanced Endoscopy Devices

Endoscopy Specialists since 1985

22134 Sherman Way
Canoga Park, CA 91303 - U.S.A.
phone: 818-227-2720 - fax: 818-227-2724
www.aed.md
E-mail: advendodev@earthlink.net