

**Advanced
Endoscopy
Devices**

AED Instrument Bridges

Patent Pending

The AED Bridge is designed to protect and elevate the tips of surgical instruments when secured onto a sterilization/instrument stringer and placed in a surgical tray.

IB100S AED Instrument Bridge - Small

Dimensions: (7" x 1" x 2")

IB100M AED Instrument Bridge - Medium


Dimensions: (11" x 1" x 2")

IB100L AED Instrument Bridge - Large

Dimensions: (16" x 1" x 2")


**Reusable &
Autoclavable**

**Manufactured from
100% Stainless Steel**



**Unique Patented Design
Protects Instrument Tips
In The Tray**

Novel bridge design helps to safely elevate the distal end of delicate or sharp surgical instrument(s) during surgery, S.P.D./C.S. processing, storage, and transportation.



- ***Unique robust design***
- ***Three sizes to fit all manufacturers stringers***
- ***Reusable and autoclavable***
- ***Environmentally conscious***
- ***Year-over-year savings by reducing repair cost***
- ***Proudly Manufactured In USA***

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md