


Advanced
Endoscopy
Devices

ENDOPRO-MDTK

3D ENDOSCOPIC VISUALIZATION SYSTEM

3D


VISUALIZATION



2D ENDOSCOPE
CAMERA SYSTEM



PRODUCT HIGHLIGHTS



1

Easy 2D/3D Mode Switch

Transition between 2D and 3D modes with a single touch.

2

Gamma Adjustment

Optimize the correlation between image contrast and brightness by adjusting the exponential relationship. Decreasing the gamma value enhances visibility in darker areas of the image.

3

3D Display Format

Choose the format that is most compatible with your 3D monitor from Line-by-line (Interlaced), Side-by-side, or Top-and-Bottom.

4

Adjustable Human Factor

Enjoy a comfortable 3D viewing experience by adjusting the 3D parameters to the preferred disparity level when discomfort or dizziness with 3D vision occurs, while maintaining accurate 3D images.

5

Tailored Algorithm for Rigid/Flexible Endoscope

Select rigid or flexible mode for optimal 3D vision tailored to your scope and procedure.

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md